

# AKTIVE DISINFECTING WIPES

**Kills Up To 99.9% of Germs,  
Viruses, and Bacteria\*  
Including SARS-CoV-2\*\***



\*KILLS: Viricidal\*, Bactericidal, Tuberculocidal, Fungicidal, Pseudomonacidal  
**AKTIVE DISINFECTING WIPES** are effective against the following pathogens in 2 minutes at 20°C/68°F, on hard, nonporous surfaces when used as directed: Escherichia coli O157: H7, Methicillin-resistant Staphylococcus aureus (MRSA), Mycobacterium tuberculosis, Pseudomonas aeruginosa, Salmonella enterica, Staphylococcus aureus, Vancomycin-resistant Enterococcus faecalis (VRE), \*Herpes Simplex Virus 2 (genital herpes virus), \*HIV-1 (AIDS virus)‡ (at an exposure time of 1 min.), \*Human Hepatitis B Virus, \*Human Hepatitis C Virus, \*Influenza A virus (A2 Japan). \*\*Kills Mycobacterium tuberculosis a harder-to-kill pathogen than SARS-CoV-2 (COVID-19).

**PREPARATION OF DISPENSER:** Remove cover from canister and discard covering seal. From center of wipes roll, pull up a wipe corner, twist into a point and thread through the tear throat in the underside of the lid cover. Pull wipe through about one (1) inch. Replace the lid cover on container, assuring a secure fit. Pull out first wipe and snap at a 90° angle. Remaining wipes feed automatically for next use. Always reseal hinged cap after each use to ensure moisture retention.

**DIRECTIONS FOR USE:** It is a violation of U.S. Federal law to use this product in a manner inconsistent with its labeling. This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. **AKTIVE DISINFECTING WIPES** may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection. This product is not for use on the inhalation therapy and/or pulmonary diagnostic equipment.

**DISINFECTION:** Thoroughly wet pre-cleaned, hard, nonporous surface with a wipe, keep wet for 2 minutes (5 minutes if fungus is suspected), and allow to air dry. Use as many wipes as needed for the treated surface to remain wet for the entire contact time. For use as a general disinfectant wipe on hard, nonporous, inanimate surfaces. You may use a wipe to pre-clean the surface to be disinfected, but thorough re-wiping with a fresh wipe is required for disinfection. Use on food-contact surfaces must be followed by a potable water rinse.

**ACTIVE INGREDIENTS:**

N-Alkyl (60% C14, 30% C16, 5% C12, 5% C18)	
Dimethyl Benzyl Ammonium Chloride.....	0.12%
N-Alkyl (68% C12, 32% C14) Dimethyl Ethyl Benzyl Ammonium Chloride.....	0.12%
Isopropyl Alcohol.....	41.58%
OTHER INGREDIENTS†.....	58.18%
<b>TOTAL:</b> .....	<b>100.00%</b>

† Does not include weight of dry wipe.

**Areas of Use:**

- \*Hotels
- \*Stadiums
- \*Schools
- \*Casinos
- \*Cities & Counties
- \*Amusement Parks
- \*Restaurants

Health	0
Flammability	2
Physical Hazard	0
Protection	X

## OMEGA INDUSTRIAL SUPPLY, INC.

1133 W 27th Street, Cheyenne, WY 82001

1-800-571-7347 [www.onlyomega.com](http://www.onlyomega.com) Fax: 307-514-0631