

MED SEAL

Drain Sealer
for Medical Environments



MED SEAL is used for the initial treatment of medical room environments such as operating, nursing or dialysis rooms. This mineral oil-based dry drain lubricant is designed to protect, lubricate and maintain dry drains in medical environments, particularly those experiencing infrequent use. With a pleasant lavender scent, **MED SEAL** serves as the first step in a comprehensive drain maintenance program. **MED SEAL** coats the interior of the pipes, preventing odors, clogs and corrosion caused by lack of water flow.

DIRECTIONS

Add sufficient water to fill the P-trap, then add 16-32 ounces of **MED SEAL** to seal in the noxious odors, reduce evaporation and deodorize. To reduce the evaporation rate even more and have the deodorizing effect last longer, do not add any water to the P-trap prior and add only **MED SEAL** to the drain. Double the below recommended maintenance dosages for these situations. Recommended dosages:

2" Drains = 6 oz

3" Drains = 32 oz

6"+ Drains = 1 gallon

| | |
|-----------------|---|
| Health | 1 |
| Flammability | 1 |
| Physical Hazard | 0 |
| Protection | B |

Areas of Use:

*Dialysis Drains

*Neptune Drains

*Floor Drains

*Hospitals

*Lab Sinks

*Medical Facilities

| | |
|-----------------------|----------------|
| Appearance..... | Purple liquid |
| Relative Density..... | 0.86 |
| Net Contents..... | 16 OZ (473 ml) |

OMEGA INDUSTRIAL SUPPLY, INC.

1133 W 27th Street, Cheyenne, WY 82001

1-800-571-7347 www.onlyomega.com Fax: 307-514-0631